



- DESCRIPTION:** COMFORT GEL A 032
- PART NUMBER:** M807 ML
- MASTER ROLL SIZE:** 37.8 MSI (THOUSAND SQ INCHES)
- GEL WIDTH:** 9 inches (22.9 cm)
- ROLL LENGTH:** 350 Linear Feet (106.68 m)
- GEL THICKNESS:** .032" +/- .004" (0.8128 mm +/-0.1016 mm)
- LINER 1:** .002" X 9.75" BLUE POLYETHYLENE (0.0508 mm X 24.76 cm)
- LINER 2:** .005" X 10" POLYESTER (0.127 mm X 25.4 cm)

MASTER ROLL




ROLL PACKAGING

- BOX SIZE & WEIGHT:** 17" X 17" X 12" (43.18cm x 43.18cm x 30.48 cm)
Approx. 75 lbs (34.1 kg)



- SHIPPING BY PALLET:** 16 boxes (Approx.54"H x 40"W x 48"D)
Approx. 1200 lbs (545.5 kg)

SLIT GEL ROLL PACKAGING:



- STORAGE:** AMBIENT CONDITIONS. RECOMMEND CONVERTING TO FINAL PACKAGING WITHIN TWELVE MONTHS OF RECEIPT TO MAINTAIN FRESHNESS AND PREVENT DRY OUT.

- ELECTRICAL:** IMPEDANCE 200 OHMS MAXIMUM AT 10 Hz
AAMI ECG PROTOCOL EC12
USING 1-INCH SQ. SILVER /SILVER CHLORIDE SUBSTRATE
VOLUME RESISTIVITY: 1000 OHM-CM MAXIMUM

- ADHESION:** PSTC- 5 TEST METHOD: 185 GRAMS/IN MINIMUM TO STAINLESS STEEL

- BIOCOMPATIBILITY:** CYTOTOXICITY STUDY USING THE ISO OVERLAY METHOD:
PER ISO 10993 NON-CYTOTOXIC, COMBINED SCORE 0.0

- ISO SKIN IRRITATION STUDY ON THE RABBIT:
NON-IRRITATING, COMBINED SCORE 0.0

- ISO SENSITIZATION STUDY ON THE GUINEA PIG:
NON-SENSITIZING, COMBINED SCORE 0.0

0127B14