

**DESCRIPTION:** COMFORT GEL B

**PART NUMBER:** M865 ML

**MASTER ROLL SIZE:** 37.8 MSI (THOUSAND SQ INCHES)

**GEL WIDTH:** 9 inches (22.9 cm)


**ROLL LENGTH:** 350 Linear Feet (106.68 m)

**GEL THICKNESS:** .032" +/- .004" (0.8128 mm +/-0.1016 mm)

**LINER 1:** .002" X 9.75" BLUE POLYETHYLENE  
(0.0508 mm X 24.76 cm)


**LINER 2:** .005" X 10" POLYESTER  
(0.127 mm X 25.4 cm)

**MASTER ROLL**




**ROLL PACKAGING**

**BOX SIZE & WEIGHT:** 17" X 17" X 12" (43.18cm x 43.18cm x 30.48 cm)  
Approx. 75 lbs (34.1 kg)



**SHIPPING BY PALLET:** 16 boxes (Approx.54"H x 40"W x 48"D)  
Approx. 1200 lbs (545.5 kg)

**SLIT GEL ROLL PACKAGING:**



**STORAGE:** AMBIENT CONDITIONS. RECOMMEND CONVERTING TO FINAL PACKAGING WITHIN TWENTY FOUR MONTHS OF RECEIPT TO MAINTAIN FRESHNESS AND PREVENT DRY OUT.

**ELECTRICAL:** IMPEDANCE 200 OHMS MAXIMUM AT 10 Hz  
AAMI ECG PROTOCOL EC12  
USING 1-INCH SQ. SILVER /SILVER CHLORIDE SUBSTRATE  
VOLUME RESISTIVITY: 1000 OHM-CM MAXIMUM

**ADHESION:** PSTC- 5 TEST METHOD: 185 GRAMS/IN MINIMUM TO STAINLESS STEEL

**BIOCOMPATIBILITY:** CYTOTOXICITY STUDY USING THE ISO OVERLAY METHOD:  
**PER ISO 10993** NON-CYTOTOXIC, COMBINED SCORE 0.0

ISO SKIN IRRITATION STUDY ON THE RABBIT:  
NON-IRRITATING, COMBINED SCORE 0.0

ISO SENSITIZATION STUDY ON THE GUINEA PIG:  
NON-SENSITIZING, COMBINED SCORE 0.0

0624C16