


DESCRIPTION: HYPERBLUE GEL B
PART NUMBER: M832B
MASTER ROLL SIZE: 32.4 MSI (THOUSAND SQ INCHES)
GEL WIDTH: 9 inches (22.9 cm)
ROLL LENGTH: 300 Linear Feet (91.44 m)
GEL THICKNESS: .045" +/- .005" (1.143 mm +/-0.127 mm)
LINER 1: .002" X 9.75" CLEAR POLYETHYLENE
(0.0508 mm X 24.76 cm)
LINER 2: .005" X 10" CLEAR POLYESTER
(0.127 mm X 25.4 cm)


MASTER ROLL




ROLL PACKAGING

BOX SIZE & WEIGHT: 17" X 17" X 12" (43.18cm x 43.18cm x 30.48 cm)
Approx. 75 lbs (34.1 kg)

SHIPPING BY PALLET: 16 boxes (Approx. 54"H x 40"W x 48"D)
Approx. 1200 lbs (545.5 kg)



SLIT GEL ROLL PACKAGING:



STORAGE: AMBIENT CONDITIONS. RECOMMEND CONVERTING TO FOIL PACKAGING WITHIN SIX MONTHS OF RECEIPT TO PREVENT DRY OUT.

ELECTRICAL: IMPEDANCE 200 OHMS MAXIMUM AT 10 Hz
AAMI ECG PROTOCOL EC12
USING 1-INCH SQ. SILVER /SILVER CHLORIDE SUBSTRATE
VOLUME RESISTIVITY: 1000 OHM-CM MAXIMUM

ADHESION: PSTC- 5 TEST METHOD: 100 GRAMS/IN MINIMUM TO STAINLESS STEEL

BIOCOMPATIBILITY: CYTOTOXICITY STUDY USING THE ISO ELUTION METHOD:
PER ISO 10993 NON-CYTOTOXIC, COMBINED SCORE 0.0
ISO SKIN IRRITATION STUDY ON THE RABBIT:
NON-IRRITATING, COMBINED SCORE 0.0
ISO SENSITIZATION STUDY ON THE GUINEA PIG:
NON-SENSITIZING, COMBINED SCORE 0.0

0127D14