



DESCRIPTION: PROMEON HI-ADHESION GEL 032
PART NUMBER: M863X
MASTER ROLL SIZE: 37.8 MSI (THOUSAND SQ INCHES)
GEL WIDTH: 9" (22.9 cm)
ROLL LENGTH: 350 Linear Feet (106.68 m)
GEL THICKNESS: .032" +/- .004"
LINER 1: .002" X 9.75" BLUE POLYETHYLENE (0.0508 mm X 24.76 cm)
LINER 2: .005" X 10" POLYESTER (0.127 mm X 25.4 cm)

MASTER ROLL




ROLL PACKAGING

BOX SIZE & WEIGHT: 17" X 17" X 12" (43.18cm x 43.18cm x 30.48 cm)
 Approx. 70 lbs (31.8 kg)



SHIPPING BY PALLET: 16 boxes (Approx. 54"H x 40"W x 48"D)
 Approx. 1120 lbs (509.1 kg)

SLIT GEL ROLL PACKAGING:



STORAGE: AMBIENT CONDITIONS. RECOMMEND CONVERTING TO FINAL PACKAGING WITHIN TWENTY FOUR MONTHS OF RECEIPT TO MAINTAIN FRESHNESS AND PREVENT DRY OUT.

ELECTRICAL: IMPEDANCE 200 OHMS MAXIMUM AT 10 Hz
 AAMI ECG PROTOCOL EC12
 USING 1-INCH SQ. SILVER /SILVER CHLORIDE SUBSTRATE
 VOLUME RESISTIVITY: @1000 Hz; 1500 OHM-CM MAXIMUM

ADHESION: PSTC- 5 TEST METHOD: 240 GRAMS/IN MINIMUM TO STAINLESS STEEL

BIOCOMPATIBILITY: CYTOTOXICITY STUDY USING THE ISO ELUTION METHOD:
PER ISO 10993 NON-CYTOTOXIC, COMBINED SCORE 0.0
 ISO SKIN IRRITATION STUDY ON THE RABBIT:
 NON-IRRITATING, COMBINED SCORE 0.0
 ISO SENSITIZATION STUDY ON THE GUINEA PIG:
 NON-SENSITIZING, COMBINED SCORE 0.0

0108B15